

# RAB9A Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01043-3

## Capture Antibody Information

**Catalog Number:**  
84129-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241372E10  
**Reactivity:**  
human  
**GenBank:**  
BC017265  
**Immunogen Catalog Number:**  
Ag1969

**Conjugate:**  
Unconjugated  
**Full name:**  
RAB9A, member RAS oncogene family  
**Gene ID:**  
9367

## Detection Antibody Information

**Catalog Number:**  
84129-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241372A6  
**Reactivity:**  
human  
**GenBank:**  
BC017265  
**Immunogen Catalog Number:**  
Ag1969

**Conjugate:**  
Unconjugated  
**Full name:**  
RAB9A, member RAS oncogene family  
**Gene ID:**  
9367

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01043-3 targets RAB9A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RAB9A Recombinant antibody, PBS Only (Capture) 84129-4-PBS (241372E10). 100 µg. Concentration 1 mg/ml.

Detection antibody: RAB9A Recombinant antibody, PBS Only (Detector) 84129-1-PBS (241372A6). 100 µg. Concentration 1 mg/ml.

Alternative RAB9A matched antibody pairs: MP01043-1, MP01043-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01043-3, RAB9A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84129-4-PBS. Detection antibody: 84129-1-PBS. Standard: Ag1969. Range: 0.313-40 ng/mL