

IL23R Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01048-1

Capture Antibody Information

Catalog Number:

84112-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240905G11

Reactivity:

human

GenBank:

NM_144701.3

Conjugate:

Unconjugated

Full name:

interleukin 23 receptor

Gene ID:

149233

Detection Antibody Information

Catalog Number:

84112-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240905D5

Reactivity:

human

GenBank:

NM_144701.3

Conjugate:

Unconjugated

Full name:

interleukin 23 receptor

Gene ID:

149233

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01048-1 targets IL23R in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL23R Recombinant antibody, PBS Only (Capture) 84112-4-PBS (240905G11). 100 µg. Concentration 1 mg/ml.

Detection antibody: IL23R Recombinant antibody, PBS Only (Detector) 84112-2-PBS (240905D5). 100 µg. Concentration 1 mg/ml.

Alternative IL23R matched antibody pairs: MP00866-1, MP01048-2, MP01048-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01048-1, IL23R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84112-4-PBS. Detection antibody: 84112-2-PBS. Standard: Eg0853. Range: 0.625-80 ng/mL