

For Research Use Only

Rat gp130/IL6ST Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01049-3

Capture Antibody Information

Catalog Number:
84113-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241340E1

Reactivity:
rat

GenBank:
NM_001008725.3

Conjugate:
Unconjugated

Full name:
interleukin 6 signal transducer

Gene ID:
25205

Detection Antibody Information

Catalog Number:
84113-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241340C9

Reactivity:
rat

GenBank:
NM_001008725.3

Conjugate:
Unconjugated

Full name:
interleukin 6 signal transducer

Gene ID:
25205

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01049-3 targets gp130/IL6ST in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Rat gp130/IL6ST Recombinant antibody, PBS Only (Capture) 84113-3-PBS (241340E1). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Rat gp130/IL6ST Recombinant antibody, PBS Only (Capture/Detector) 84113-2-PBS (241340C9). 100 μ g. Concentration 1 mg/ml.

Alternative gp130/IL6ST matched antibody pairs: MP00324-1, MP00324-2, MP00324-3, MP01049-1, MP01049-2, MP50086-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptgcn.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01049-3, RAT gp130 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84113-3-PBS. Detection antibody: 84113-2-PBS. Standard: Eg0857. Range: 0.625-80 ng/mL.