

For Research Use Only

CD155/PVR Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01053-1

Capture Antibody Information

Catalog Number:
84138-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241405B2

Reactivity:
human

GenBank:
BC015542

Conjugate:
Unconjugated

Full name:
poliovirus receptor

Gene ID:
5817

Detection Antibody Information

Catalog Number:
84138-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241405D2

Reactivity:
human

GenBank:
BC015542

Conjugate:
Unconjugated

Full name:
poliovirus receptor

Gene ID:
5817

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01053-1 targets CD155/PVR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD155/PVR Recombinant antibody, PBS Only (Capture) 84138-1-PBS (241405B2). 100 μ g. Concentration 1 mg/ml.

Detection antibody: CD155/PVR Recombinant antibody, PBS Only (Detector) 84138-2-PBS (241405D2). 100 μ g. Concentration 1 mg/ml.

Alternative CD155/PVR matched antibody pairs: MP00705-1, MP00705-2, MP00705-3, MP01053-2, MP01053-3, MP01053-4, MP50402-1, MP50402-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01053-1, CD155/PVR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84138-1-PBS. Detection antibody: 84138-2-PBS. Standard: Eg0971. Range: 0.313-40 ng/mL.