

For Research Use Only

CCDC134 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01057-1

Capture Antibody Information

Catalog Number:
84106-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241212C8
Reactivity:
human
GenBank:
BC017693
Immunogen Catalog Number:
Ag21990

Conjugate:
Unconjugated
Full name:
coiled-coil domain containing 134
Gene ID:
79879

Detection Antibody Information

Catalog Number:
84106-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241212F11
Reactivity:
human
GenBank:
BC017693
Immunogen Catalog Number:
Ag21990

Conjugate:
Unconjugated
Full name:
coiled-coil domain containing 134
Gene ID:
79879

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01057-1 targets CCDC134 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: CCDC134 Recombinant antibody, PBS Only (Capture/Detector) 84106-1-PBS (241212C8). 100 μ g. Concentration 1 mg/ml.
Detection antibody: CCDC134 Recombinant antibody, PBS Only (Detector) 84106-3-PBS (241212F11). 100 μ g. Concentration 1 mg/ml.
Alternative CCDC134 matched antibody pairs: MP01057-2, MP01057-3
Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01057-1, CCDC134 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84106-1-PBS. Detection antibody: 84106-3-PBS. Standard: Ag21990. Range: 0.625-80 ng/mL.