

ABCC10 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01067-2

Capture Antibody Information

Catalog Number:
84164-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241429F2
Reactivity:
human
GenBank:
NM_001198934.1
Immunogen Catalog Number:
Ag33719

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family C (CFTR/MRP), member 10
Gene ID:
89845

Detection Antibody Information

Catalog Number:
84164-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241429E6
Reactivity:
human
GenBank:
NM_001198934.1
Immunogen Catalog Number:
Ag33719

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family C (CFTR/MRP), member 10
Gene ID:
89845

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01067-2 targets ABCC10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ABCC10 Recombinant antibody, PBS Only (Capture/Detector) 84164-3-PBS (241429F2). 100 µg. Concentration 1 mg/ml.

Detection antibody: ABCC10 Recombinant antibody, PBS Only (Detector) 84164-2-PBS (241429E6). 100 µg. Concentration 1 mg/ml.

Alternative ABCC10 matched antibody pairs: MP01067-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01067-2, ABCC10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84164-3-PBS. Detection antibody: 84164-2-PBS. Standard: Ag33719. Range: 0.781-100 ng/mL