

For Research Use Only

# PCSK9 Recombinant Matched Antibody Pair, PBS Only



Catalog Number: MP01071-3

## Capture Antibody Information

Catalog Number:  
84172-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241395G1

Reactivity:  
human

GenBank:  
NM\_174936.4

Conjugate:  
Unconjugated

Full name:  
proprotein convertase  
subtilisin/kexin type 9

Gene ID:  
255738

## Detection Antibody Information

Catalog Number:  
84172-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241395E11

Reactivity:  
human

GenBank:  
NM\_174936.4

Conjugate:  
Unconjugated

Full name:  
proprotein convertase  
subtilisin/kexin type 9

Gene ID:  
255738

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01071-3 targets PCSK9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCSK9 Recombinant antibody, PBS Only (Capture) 84172-4-PBS (241395G1). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: PCSK9 Recombinant antibody, PBS Only (Detector) 84172-3-PBS (241395E11). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative PCSK9 matched antibody pairs: MP00059-1, MP00059-2, MP00059-3, MP01071-1, MP01071-2, MP50055-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01071-3, PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84172-4-PBS. Detection antibody: 84172-3-PBS. Standard: Eg0362. Range: 1.563-100 ng/mL.