

For Research Use Only

# Nepriylsin/CD10 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01076-3

## Capture Antibody Information

Catalog Number:

84181-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241323A4

Reactivity:

human

GenBank:

BC101658

Conjugate:

Unconjugated

Full name:

membrane metallo-endopeptidase

Gene ID:

4311

## Detection Antibody Information

Catalog Number:

84181-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241322A8

Reactivity:

human

GenBank:

BC101658

Conjugate:

Unconjugated

Full name:

membrane metallo-endopeptidase

Gene ID:

4311

## Applications

Tested Applications:

Cytometric bead array

Range:

0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01076-3 targets Nepriylsin/CD10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Nepriylsin/CD10 Recombinant antibody, PBS Only (Capture/Detector) 84181-2-PBS (241323A4). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: Nepriylsin/CD10 Recombinant antibody, PBS Only (Detector) 84181-1-PBS (241322A8). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative Nepriylsin/CD10 matched antibody pairs: MP01076-1, MP01076-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01076-3, MME,CD10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84181-2-PBS. Detection antibody: 84181-1-PBS. Standard: RP02332. Range: 0.781-100 ng/mL