

# ERBB2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01080-2

## Capture Antibody Information

**Catalog Number:**  
84182-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241409F9  
**Reactivity:**  
human  
**GenBank:**  
BC156755

**Conjugate:**  
Unconjugated  
**Full name:**  
v-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian)  
**Gene ID:**  
2064

## Detection Antibody Information

**Catalog Number:**  
84182-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241409D2  
**Reactivity:**  
human  
**GenBank:**  
BC156755

**Conjugate:**  
Unconjugated  
**Full name:**  
v-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian)  
**Gene ID:**  
2064

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01080-2 targets ERBB2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ERBB2 Recombinant antibody, PBS Only (Capture) 84182-2-PBS (241409F9). 100 µg. Concentration 1 mg/mL.

Detection antibody: ERBB2 Recombinant antibody, PBS Only (Detector) 84182-1-PBS (241409D2). 100 µg. Concentration 1 mg/mL.

Alternative ERBB2 matched antibody pairs: MP01080-1, MP01080-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01080-2, ERBB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84182-2-PBS. Detection antibody: 84182-1-PBS. Standard: RP02350. Range: 0.625-80 ng/mL