

IL-7 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01089-2

Capture Antibody Information

Catalog Number:
84165-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241000E1

Reactivity:
human

GenBank:
NM_000880.4

Conjugate:
Unconjugated

Full name:
interleukin 7

Gene ID:
3574

Detection Antibody Information

Catalog Number:
84165-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241000E2

Reactivity:
human

GenBank:
NM_000880.4

Conjugate:
Unconjugated

Full name:
interleukin 7

Gene ID:
3574

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01089-2 targets IL-7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL-7 Recombinant antibody, PBS Only (Capture) 84165-1-PBS (241000E1). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL-7 Recombinant antibody, PBS Only (Detector) 84165-2-PBS (241000E2). 100 µg. Concentration 1 mg/mL.

Alternative IL-7 matched antibody pairs: MP01089-1, MP01089-3, MP50135-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01089-2, IL-7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84165-1-PBS. Detection antibody: 84165-2-PBS. Standard: Eg0830. Range: 0.313-20 ng/mL