

# PDCD4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01098-1

## Capture Antibody Information

**Catalog Number:**  
84162-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241463E3  
**Reactivity:**  
human  
**GenBank:**  
BC026104  
**Immunogen Catalog Number:**  
Ag3288

**Conjugate:**  
Unconjugated  
**Full name:**  
programmed cell death 4 (neoplastic transformation inhibitor)  
**Gene ID:**  
27250

## Detection Antibody Information

**Catalog Number:**  
84162-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241463A8  
**Reactivity:**  
human  
**GenBank:**  
BC026104  
**Immunogen Catalog Number:**  
Ag3288

**Conjugate:**  
Unconjugated  
**Full name:**  
programmed cell death 4 (neoplastic transformation inhibitor)  
**Gene ID:**  
27250

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01098-1 targets PDCD4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDCD4 Recombinant antibody, PBS Only (Capture) 84162-3-PBS (241463E3). 100 µg. Concentration 1 mg/ml.

Detection antibody: PDCD4 Recombinant antibody, PBS Only (Capture/Detector) 84162-1-PBS (241463A8). 100 µg. Concentration 1 mg/ml.

Alternative PDCD4 matched antibody pairs: MP01098-2, MP01098-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01098-1, PDCD4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84162-3-PBS. Detection antibody: 84162-1-PBS. Standard: Ag3288. Range: 0.781-100 ng/mL.