

For Research Use Only

ANO1,DOG1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01115-1

Capture Antibody Information

Catalog Number:
83567-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241464H5
Reactivity:
human
GenBank:
BC027590
Immunogen Catalog Number:
Ag3320

Conjugate:
Unconjugated
Full name:
anoctamin 1, calcium activated chloride channel
Gene ID:
55107

Detection Antibody Information

Catalog Number:
83567-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241464C11
Reactivity:
human
GenBank:
BC027590
Immunogen Catalog Number:
Ag3320

Conjugate:
Unconjugated
Full name:
anoctamin 1, calcium activated chloride channel
Gene ID:
55107

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01115-1 targets ANO1,DOG1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ANO1,DOG1 Recombinant antibody, PBS Only (Capture) 83567-5-PBS (241464H5). 100 µg. Concentration 1 mg/ml.

Detection antibody: ANO1,DOG1 Recombinant antibody, PBS Only (Detector) 83567-3-PBS (241464C11). 100 µg. Concentration 1 mg/ml.

Alternative ANO1,DOG1 matched antibody pairs: MP01115-2, MP50123-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01115-1, ANO1,DOG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83567-5-PBS. Detection antibody: 83567-3-PBS. Standard: Ag3320. Range: 0.313-40 ng/mL.