

For Research Use Only

CD122/IL-2RB Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01123-1

Capture Antibody Information

Catalog Number: 84229-3-PBS	Clone ID: 241408B1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: interleukin 2 receptor, beta
Isotype: IgG	GenBank: BC025691	Gene ID: 3560
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3929	

Detection Antibody Information

Catalog Number: 84229-2-PBS	Clone ID: 241408A11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: interleukin 2 receptor, beta
Isotype: IgG	GenBank: BC025691	Gene ID: 3560
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3929	

Applications

Tested Applications: Cytometric bead array	Range: 0156-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01123-1 targets CD122/IL-2RB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD122/IL-2RB Recombinant antibody, PBS Only (Capture) 84229-3-PBS (241408B1). 100 µg. Concentration 1 mg/ml.

Detection antibody: CD122/IL-2RB Recombinant antibody, PBS Only (Detector) 84229-2-PBS (241408A11). 100 µg. Concentration 1 mg/ml.

Alternative CD122/IL-2RB matched antibody pairs: MP01123-2, MP01123-3, MP50451-1, MP50451-2, MP50451-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01123-1, CD122/IL2RB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84229-3-PBS. Detection antibody: 84229-2-PBS. Standard: Ag3929. Range: 0.156-20 ng/mL.