

FGFR2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01128-3

Capture Antibody Information

Catalog Number:
84205-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241383B5

Reactivity:
human

GenBank:
NM_000141.4

Conjugate:
Unconjugated

Full name:
fibroblast growth factor receptor 2

Gene ID:
2263

Detection Antibody Information

Catalog Number:
84205-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241383G9

Reactivity:
human

GenBank:
NM_000141.4

Conjugate:
Unconjugated

Full name:
fibroblast growth factor receptor 2

Gene ID:
2263

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01128-3 targets FGFR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FGFR2 Recombinant antibody, PBS Only (Capture/Detector) 84205-2-PBS (241383B5). 100 µg. Concentration 1 mg/ml.

Detection antibody: FGFR2 Recombinant antibody, PBS Only (Detector) 84205-4-PBS (241383G9). 100 µg. Concentration 1 mg/ml.

Alternative FGFR2 matched antibody pairs: MP01128-1, MP01128-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01128-3, FGFR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84205-2-PBS. Detection antibody: 84205-4-PBS. Standard: Eg1265. Range: 0.625-80 ng/mL