

For Research Use Only

LIF Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01130-1

Capture Antibody Information

Catalog Number:
84232-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241519A7
Reactivity:
human
GenBank:
BC093733

Conjugate:
Unconjugated
Full name:
leukemia inhibitory factor
(cholinergic differentiation factor)
Gene ID:
3976

Detection Antibody Information

Catalog Number:
84232-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241519A6
Reactivity:
human
GenBank:
BC093733

Conjugate:
Unconjugated
Full name:
leukemia inhibitory factor
(cholinergic differentiation factor)
Gene ID:
3976

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01130-1 targets LIF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LIF Recombinant antibody, PBS Only (Capture) 84232-3-PBS (241519A7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: LIF Recombinant antibody, PBS Only (Detector) 84232-2-PBS (241519A6). 100 μ g. Concentration 1 mg/mL.

Alternative LIF matched antibody pairs: MP01130-2, MP01130-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01130-1, LIF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84232-3-PBS. Detection antibody: 84232-2-PBS. Standard: Eg0111. Range: 0.625-80 ng/mL