

# KLHDC3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01151-2

## Capture Antibody Information

<b>Catalog Number:</b> 84217-4-PBS	<b>Clone ID:</b> 241475G5	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> kelch domain containing 3
<b>Isotype:</b> IgG	<b>GenBank:</b> BC041793	<b>Gene ID:</b> 116138
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag4879	

## Detection Antibody Information

<b>Catalog Number:</b> 84217-2-PBS	<b>Clone ID:</b> 241475B12	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> kelch domain containing 3
<b>Isotype:</b> IgG	<b>GenBank:</b> BC041793	<b>Gene ID:</b> 116138
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag4879	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.781-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP01151-2 targets KLHDC3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KLHDC3 Recombinant antibody, PBS Only (Capture) 84217-4-PBS (241475G5). 100 µg. Concentration 1 mg/ml.

Detection antibody: KLHDC3 Recombinant antibody, PBS Only (Capture/Detector) 84217-2-PBS (241475B12). 100 µg. Concentration 1 mg/ml.

Alternative KLHDC3 matched antibody pairs: MP01151-1, MP01151-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01151-2, KLHDC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84217-4-PBS. Detection antibody: 84217-2-PBS. Standard: Ag4879. Range: 0.781-100 ng/mL.