

TNFRSF13B/CD267 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01167-1

Capture Antibody Information

Catalog Number:
84266-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241317B5
Reactivity:
human
GenBank:
BC109392

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 13B
Gene ID:
23495

Detection Antibody Information

Catalog Number:
98100-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241317D8
Reactivity:
human
GenBank:
BC109392

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 13B
Gene ID:
23495

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01167-1 targets TNFRSF13B/CD267 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFRSF13B Recombinant antibody, PBS Only (Capture) 84266-1-PBS (241317B5). 100 µg. Concentration 1 mg/ml.

Detection antibody: Anti-Human TNFRSF13B/CD267 Rabbit Recombinant Antibody, PBS Only 98100-1-PBS (241317D8). 100 µg. Concentration 1 mg/ml.

Alternative TNFRSF13B/CD267 matched antibody pairs: MP01167-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01167-1, TNFRSF13B/CD267 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84266-1-PBS. Detection antibody: 98100-1-PBS. Standard: RP02325. Range: 0.313-40 ng/mL.