

BTLA/CD272 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01168-1

Capture Antibody Information

Catalog Number:
84262-5-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241319B1

Reactivity:
human

GenBank:
BC107091

Conjugate:
Unconjugated

Full name:
B and T lymphocyte associated

Gene ID:
151888

Detection Antibody Information

Catalog Number:
84262-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241319C12

Reactivity:
human

GenBank:
BC107091

Conjugate:
Unconjugated

Full name:
B and T lymphocyte associated

Gene ID:
151888

Applications

Tested Applications:
Cytometric bead array

Range:
0.078-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01168-1 targets BTLA/CD272 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BTLA/CD272 Recombinant antibody, PBS Only (Capture) 84262-5-PBS (241319B1). 100 μ g. Concentration 1 mg/ml.

Detection antibody: BTLA/CD272 Recombinant antibody, PBS Only (Detector) 84262-4-PBS (241319C12). 100 μ g. Concentration 1 mg/ml.

Alternative BTLA/CD272 matched antibody pairs: MP01168-2, MP01168-3, MP01168-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01168-1, BTLA/CD272 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84262-5-PBS. Detection antibody: 84262-4-PBS. Standard: RP02326. Range: 0.078-10 ng/mL.