

BTLA/CD272 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01168-2

Capture Antibody Information

Catalog Number:

98096-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241318H5

Reactivity:

human

GenBank:

BC107091

Conjugate:

Unconjugated

Full name:

B and T lymphocyte associated

Gene ID:

151888

Detection Antibody Information

Catalog Number:

84262-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241319C10

Reactivity:

human

GenBank:

BC107091

Conjugate:

Unconjugated

Full name:

B and T lymphocyte associated

Gene ID:

151888

Applications

Tested Applications:

Cytometric bead array

Range:

0.078-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01168-2 targets BTLA/CD272 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Anti-Human BTLA/CD272 Rabbit Recombinant Antibody, PBS Only 98096-1-PBS (241318H5). 100 μ g. Concentration 1 mg/ml.

Detection antibody: BTLA/CD272 Recombinant antibody, PBS Only (Detector) 84262-2-PBS (241319C10). 100 μ g. Concentration 1 mg/ml.

Alternative BTLA/CD272 matched antibody pairs: MP01168-1, MP01168-3, MP01168-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01168-2, BTLA/CD272 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 98096-1-PBS. Detection antibody: 84262-2-PBS. Standard: RP02326. Range: 0.078-10 ng/mL.