

CD72 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01170-1

Capture Antibody Information

Catalog Number:
84270-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241301D3

Reactivity:
human

GenBank:
BC030227

Conjugate:
Unconjugated

Full name:
CD72 molecule

Gene ID:
971

Detection Antibody Information

Catalog Number:
84270-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241300G3

Reactivity:
human

GenBank:
BC030227

Conjugate:
Unconjugated

Full name:
CD72 molecule

Gene ID:
971

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01170-1 targets CD72 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD72 Recombinant antibody, PBS Only (Capture) 84270-2-PBS (241301D3). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD72 Recombinant antibody, PBS Only (Capture/Detector) 84270-4-PBS (241300G3). 100 µg. Concentration 1 mg/mL.

Alternative CD72 matched antibody pairs: MP01170-2, MP01170-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01170-1, CD72 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84270-2-PBS. Detection antibody: 84270-4-PBS. Standard: RP02318. Range: 0.313-40 ng/mL