

For Research Use Only



Noggin Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01174-1

Capture Antibody Information

Catalog Number:
84283-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241511C9

Reactivity:
human

GenBank:
BC034027

Immunogen Catalog Number:
Ag32755

Conjugate:
Unconjugated

Full name:
noggin

Gene ID:
9241

Detection Antibody Information

Catalog Number:
84283-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241511A2

Reactivity:
human

GenBank:
BC034027

Immunogen Catalog Number:
Ag32755

Conjugate:
Unconjugated

Full name:
noggin

Gene ID:
9241

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01174-1 targets Noggin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NOG Recombinant antibody, PBS Only (Capture) 84283-1-PBS (241511C9). 100 μ g. Concentration 1 mg/mL.

Detection antibody: NOG Recombinant antibody, PBS Only (Capture/Detector) 84283-2-PBS (241511A2). 100 μ g. Concentration 1 mg/mL.

Alternative Noggin matched antibody pairs: MP01174-2, MP01174-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01174-1, Noggin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84283-1-PBS. Detection antibody: 84283-2-PBS. Standard: Ag32755. Range: 0.313-40 ng/mL