

# TOR1B Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01175-1

## Capture Antibody Information

Catalog Number:  
84267-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241574F4

Reactivity:  
human

GenBank:  
BC015578

Immunogen Catalog Number:  
Ag16641

Conjugate:  
Unconjugated

Full name:  
torsin family 1, member B (torsin B)

Gene ID:  
27348

## Detection Antibody Information

Catalog Number:  
84267-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241574A7

Reactivity:  
human

GenBank:  
BC015578

Immunogen Catalog Number:  
Ag16641

Conjugate:  
Unconjugated

Full name:  
torsin family 1, member B (torsin B)

Gene ID:  
27348

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01175-1 targets TOR1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TOR1B Recombinant antibody, PBS Only (Capture) 84267-1-PBS (241574F4). 100 µg. Concentration 1 mg/ml.

Detection antibody: TOR1B Recombinant antibody, PBS Only (Detector) 84267-3-PBS (241574A7). 100 µg. Concentration 1 mg/ml.

Alternative TOR1B matched antibody pairs: MP01175-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:


Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01175-1, TOR1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84267-1-PBS. Detection antibody: 84267-3-PBS. Standard: Ag16641. Range: 0.313-40 ng/mL.