

TOR1B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01175-2

Capture Antibody Information

Catalog Number:
84267-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241574G4
Reactivity:
human
GenBank:
BC015578
Immunogen Catalog Number:
Ag16641

Conjugate:
Unconjugated
Full name:
torsin family 1, member B (torsin B)
Gene ID:
27348

Detection Antibody Information

Catalog Number:
84267-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241574A7
Reactivity:
human
GenBank:
BC015578
Immunogen Catalog Number:
Ag16641

Conjugate:
Unconjugated
Full name:
torsin family 1, member B (torsin B)
Gene ID:
27348

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01175-2 targets TOR1B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TOR1B Recombinant antibody, PBS Only (Capture) 84267-2-PBS (241574G4). 100 µg. Concentration 1 mg/ml.

Detection antibody: TOR1B Recombinant antibody, PBS Only (Detector) 84267-3-PBS (241574A7). 100 µg. Concentration 1 mg/ml.

Alternative TOR1B matched antibody pairs: MP01175-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01175-2, TOR1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84267-2-PBS. Detection antibody: 84267-3-PBS. Standard: Ag16641. Range: 0.313-40 ng/mL.