

Complement component C7/C7 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01179-1

Capture Antibody Information

Catalog Number:
84263-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241551E11
Reactivity:
human
GenBank:
NM_000587.4

Conjugate:
Unconjugated
Full name:
complement component 7
Gene ID:
730

Detection Antibody Information

Catalog Number:
84263-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241551D6
Reactivity:
human
GenBank:
NM_000587.4

Conjugate:
Unconjugated
Full name:
complement component 7
Gene ID:
730

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01179-1 targets Complement component C7/C7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Complement component C7/C7 Recombinant antibody, PBS Only (Capture/Detector) 84263-1-PBS (241551E11). 100 µg. Concentration 1 mg/ml.

Detection antibody: Complement component C7/C7 Recombinant antibody, PBS Only (Detector) 84263-5-PBS (241551D6). 100 µg. Concentration 1 mg/ml.

Alternative Complement component C7/C7 matched antibody pairs: MP01179-2, MP01179-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01179-1, Complement component C7/C7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84263-1-PBS. Detection antibody: 84263-5-PBS. Standard: Eg1202. Range: 0.625-80 ng/mL