

# ATPBD3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP01189-1

## Capture Antibody Information

Catalog Number:  
84318-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241630E9

Reactivity:  
human

GenBank:  
NM\_145232

Immunogen Catalog Number:  
Ag35037

Conjugate:  
Unconjugated

Full name:  
ATP binding domain 3

Gene ID:  
90353

## Detection Antibody Information

Catalog Number:  
84318-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241630G1

Reactivity:  
human

GenBank:  
NM\_145232

Immunogen Catalog Number:  
Ag35037

Conjugate:  
Unconjugated

Full name:  
ATP binding domain 3

Gene ID:  
90353

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01189-1 targets ATPBD3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ATPBD3 Recombinant antibody, PBS Only (Capture) 84318-2-PBS (241630E9). 100 µg. Concentration 1 mg/ml.

Detection antibody: ATPBD3 Recombinant antibody, PBS Only (Capture/Detector) 84318-3-PBS (241630G1). 100 µg. Concentration 1 mg/ml.

Alternative ATPBD3 matched antibody pairs: MP01189-2, MP01189-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01189-1, ATPBD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84318-2-PBS. Detection antibody: 84318-3-PBS. Standard: Ag35037 Range: 0.625-80 ng/mL.