

For Research Use Only



RINT1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01190-2

Capture Antibody Information

Catalog Number: 84290-4-PBS	Clone ID: 241592D4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: RAD50 interactor 1
Isotype: IgG	GenBank: BC068483	Gene ID: 60561
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6145	

Detection Antibody Information

Catalog Number: 84290-5-PBS	Clone ID: 241592A8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: RAD50 interactor 1
Isotype: IgG	GenBank: BC068483	Gene ID: 60561
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6145	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP01190-2 targets RINT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RINT1 Recombinant antibody, PBS Only (Capture) 84290-4-PBS (241592D4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: RINT1 Recombinant antibody, PBS Only (Detector) 84290-5-PBS (241592A8). 100 μ g. Concentration 1 mg/ml.

Alternative RINT1 matched antibody pairs: MP01190-1, MP01190-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01190-2, RINT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84290-4-PBS. Detection antibody: 84290-5-PBS. Standard: Ag6145. Range: 0.625-80 ng/mL