

For Research Use Only

TMEM98 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01195-2

Capture Antibody Information

| | | |
|---|--|---|
| Catalog Number: 84284-1-PBS | Clone ID: 241581B9 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: transmembrane protein 98 |
| Isotype: IgG | GenBank: BC000526 | Gene ID: 26022 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag6492 | |

Detection Antibody Information

| | | |
|---|--|---|
| Catalog Number: 84284-3-PBS | Clone ID: 241581A2 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: transmembrane protein 98 |
| Isotype: IgG | GenBank: BC000526 | Gene ID: 26022 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag6492 | |

Applications

| | | |
|--|---|---|
| Tested Applications: Cytometric bead array | Range: 0.313-40 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|--|---|---|

Product Information

MP01195-2 targets TMEM98 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMEM98 Recombinant antibody, PBS Only (Capture) 84284-1-PBS (241581B9). 100 µg. Concentration 1 mg/ml.

Detection antibody: TMEM98 Recombinant antibody, PBS Only (Detector) 84284-3-PBS (241581A2). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01195-2, TMEM98 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84284-1-PBS. Detection antibody: 84284-3-PBS. Standard: Ag6492. Range:0.313-40 ng/mL.