

# NME3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01200-2

## Capture Antibody Information

**Catalog Number:**  
84278-5-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241585H7  
**Reactivity:**  
human  
**GenBank:**  
BC000250  
**Immunogen Catalog Number:**  
Ag7284

**Conjugate:**  
Unconjugated  
**Full name:**  
non-metastatic cells 3, protein expressed in  
**Gene ID:**  
4832

## Detection Antibody Information

**Catalog Number:**  
84278-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241585D6  
**Reactivity:**  
human  
**GenBank:**  
BC000250  
**Immunogen Catalog Number:**  
Ag7284

**Conjugate:**  
Unconjugated  
**Full name:**  
non-metastatic cells 3, protein expressed in  
**Gene ID:**  
4832

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01200-2 targets NME3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NME3 Recombinant antibody, PBS Only (Capture) 84278-5-PBS (241585H7). 100 µg. Concentration 1 mg/mL.

Detection antibody: NME3 Recombinant antibody, PBS Only (Detector) 84278-1-PBS (241585D6). 100 µg. Concentration 1 mg/mL.

Alternative NME3 matched antibody pairs: MP01200-1, MP01200-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01200-2, NME3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84278-5-PBS. Detection antibody: 84278-1-PBS. Standard: Ag7284. Range: 0.313-40 ng/mL