

For Research Use Only



ALKBH3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01211-2

Capture Antibody Information

Catalog Number:
84314-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241650F10
Reactivity:
human, mouse, rat
GenBank:
BC015155
Immunogen Catalog Number:
Ag2982

Conjugate:
Unconjugated
Full name:
alkB, alkylation repair homolog 3 (E. coli)
Gene ID:
221120

Detection Antibody Information

Catalog Number:
84314-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241650D3
Reactivity:
human
GenBank:
BC015155
Immunogen Catalog Number:
Ag2982

Conjugate:
Unconjugated
Full name:
alkB, alkylation repair homolog 3 (E. coli)
Gene ID:
221120

Applications

Tested Applications:
Sandwich ELISA

Range:
0.156-10 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01211-2 targets ALKBH3 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: ALKBH3 Recombinant antibody, PBS Only (Capture) 84314-2-PBS (241650F10). 100 µg. Concentration 1 mg/ml.

Detection antibody: ALKBH3 Recombinant antibody, PBS Only (Detector) 84314-4-PBS (241650D3). 100 µg. Concentration 1 mg/ml.

Alternative ALKBH3 matched antibody pairs: MP01211-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01211-2, Human ALKBH3 Recombinant Matched Antibody Pair - PBS only. 84314-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2982. 84314-4-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.



The mean ALKBH3 concentration was determined to be 6.80 ng/mL in Jurkat cell extract based on a 1.50 mg/mL extract load.