

ACAP1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01212-1

Capture Antibody Information

Catalog Number:
82967-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230370G6
Reactivity:
human
GenBank:
BC018543
Immunogen Catalog Number:
Ag34309

Conjugate:
Unconjugated
Full name:
ArfGAP with coiled-coil, ankyrin repeat and PH domains 1
Gene ID:
9744

Detection Antibody Information

Catalog Number:
82967-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230370G9
Reactivity:
human
GenBank:
BC018543
Immunogen Catalog Number:
Ag34309

Conjugate:
Unconjugated
Full name:
ArfGAP with coiled-coil, ankyrin repeat and PH domains 1
Gene ID:
9744

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01212-1 targets ACAP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ACAP1 Recombinant antibody, PBS Only (Capture) 82967-2-PBS (230370G6). 100 µg. Concentration 1 mg/ml.

Detection antibody: ACAP1 Recombinant antibody, PBS Only (Detector) 82967-3-PBS (230370G9). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01212-1, ACAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82967-2-PBS. Detection antibody: 82967-3-PBS. Standard: Ag34309. Range: 1.25-20 ng/mL