

FTSJ3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01213-2

Capture Antibody Information

Catalog Number: 84311-4-PBS	Clone ID: 241507F6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: FtsJ homolog 3 (E. coli)
Isotype: IgG	GenBank: BC036710	Gene ID: 117246
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9835	

Detection Antibody Information

Catalog Number: 84311-2-PBS	Clone ID: 241507D1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: FtsJ homolog 3 (E. coli)
Isotype: IgG	GenBank: BC036710	Gene ID: 117246
Purification Method: Protein A purification	Immunogen Catalog Number: Ag9835	

Applications

Tested Applications: Cytometric bead array	Range: 0.156-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	---------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP01213-2 targets FTSJ3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FTSJ3 Recombinant antibody, PBS Only (Capture) 84311-4-PBS (241507F6). 100 µg. Concentration 1 mg/ml.

Detection antibody: FTSJ3 Recombinant antibody, PBS Only (Detector) 84311-2-PBS (241507D1). 100 µg. Concentration 1 mg/ml.

Alternative FTSJ3 matched antibody pairs: MP01213-1, MP01213-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01213-2, FTSJ3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84311-4-PBS. Detection antibody: 84311-2-PBS. Standard: Ag9835. Range: 0.156-20 ng/mL