

# FTSJ3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01213-3

## Capture Antibody Information

<b>Catalog Number:</b> 84311-3-PBS	<b>Clone ID:</b> 241507E3	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> FtsJ homolog 3 (E. coli)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC036710	<b>Gene ID:</b> 117246
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag9835	

## Detection Antibody Information

<b>Catalog Number:</b> 84311-1-PBS	<b>Clone ID:</b> 241507C3	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> FtsJ homolog 3 (E. coli)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC036710	<b>Gene ID:</b> 117246
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag9835	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.156-20 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP01213-3 targets FTSJ3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FTSJ3 Recombinant antibody, PBS Only (Capture) 84311-3-PBS (241507E3). 100 µg. Concentration 1 mg/ml.

Detection antibody: FTSJ3 Recombinant antibody, PBS Only (Detector) 84311-1-PBS (241507C3). 100 µg. Concentration 1 mg/ml.

Alternative FTSJ3 matched antibody pairs: MP01213-1, MP01213-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP01213-3, FTSJ3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84311-3-PBS. Detection antibody: 84311-1-PBS. Standard: Ag9835. Range: 0.156-20 ng/mL