

For Research Use Only

Alpha-2-Macroglobulin Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01224-1

Capture Antibody Information

Catalog Number:
84347-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241649D11

Reactivity:
human

GenBank:
NM_000014.6

Conjugate:
Unconjugated

Full name:
alpha-2-macroglobulin

Gene ID:
2

Detection Antibody Information

Catalog Number:
84347-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241649F1

Reactivity:
human

GenBank:
NM_000014.6

Conjugate:
Unconjugated

Full name:
alpha-2-macroglobulin

Gene ID:
2

Applications

Tested Applications:
Cytometric bead array

Range:
0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01224-1 targets Alpha-2-Macroglobulin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Alpha-2-Macroglobulin Recombinant antibody, PBS Only (Capture) 84347-1-PBS (241649D11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: Alpha-2-Macroglobulin Recombinant antibody, PBS Only (Detector) 84347-2-PBS (241649F1). 100 μ g. Concentration 1 mg/mL.

Alternative Alpha-2-Macroglobulin matched antibody pairs: MP01224-2, MP01224-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01224-1, Alpha-2-macroglobulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84347-1-PBS. Detection antibody: 84347-2-PBS. Standard:null. Range: null.