

For Research Use Only

POLR3C Recombinant Matched Antibody Pair, PBS Only



Catalog Number: MP01234-3

Capture Antibody Information

Catalog Number:
84367-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241402A4

Reactivity:
human

GenBank:
BC002586

Immunogen Catalog Number:
Ag27311

Conjugate:
Unconjugated

Full name:
polymerase (RNA) III (DNA directed)
polypeptide C (62kD)

Gene ID:
10623

Detection Antibody Information

Catalog Number:
84367-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241402C9

Reactivity:
human

GenBank:
BC002586

Immunogen Catalog Number:
Ag27311

Conjugate:
Unconjugated

Full name:
polymerase (RNA) III (DNA directed)
polypeptide C (62kD)

Gene ID:
10623

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01234-3 targets POLR3C in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: POLR3C Recombinant antibody, PBS Only (Capture) 84367-1-PBS (241402A4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: POLR3C Recombinant antibody, PBS Only (Detector) 84367-3-PBS (241402C9). 100 μ g. Concentration 1 mg/ml.

Alternative POLR3C matched antibody pairs: MP01234-1, MP01234-2, MP01234-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01234-3, POLR3C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84367-1-PBS. Detection antibody: 84367-3-PBS. Standard: Ag27311. Range: 0.781-100 ng/mL