

For Research Use Only



# PDGFRA Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01256-1

## Capture Antibody Information

**Catalog Number:**  
84383-5-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241659E9  
**Reactivity:**  
human  
**GenBank:**  
BC015186

**Conjugate:**  
Unconjugated  
**Full name:**  
platelet-derived growth factor receptor, alpha polypeptide  
**Gene ID:**  
5156

## Detection Antibody Information

**Catalog Number:**  
98141-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241658H4  
**Reactivity:**  
human  
**GenBank:**  
BC015186

**Conjugate:**  
Unconjugated  
**Full name:**  
platelet-derived growth factor receptor, alpha polypeptide  
**Gene ID:**  
5156

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01256-1 targets PDGFRA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDGFRA Recombinant antibody, PBS Only (Capture) 84383-5-PBS (241659E9). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: Anti-Human PDGFRA Rabbit Recombinant Antibody, PBS Only 98141-1-PBS (241658H4). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative PDGFRA matched antibody pairs: MP01256-2, MP01256-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01256-1, PDGFRA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84383-5-PBS. Detection antibody: 98141-1-PBS. Standard: RP00003. Range: 0.781-100 ng/mL