

For Research Use Only

CD89 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01257-1

Capture Antibody Information

Catalog Number:
84384-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241426D1

Reactivity:
human

GenBank:
BC027953

Conjugate:
Unconjugated

Full name:
Fc fragment of IgA, receptor for

Gene ID:
2204

Detection Antibody Information

Catalog Number:
84384-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241425F6

Reactivity:
human

GenBank:
BC027953

Conjugate:
Unconjugated

Full name:
Fc fragment of IgA, receptor for

Gene ID:
2204

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01257-1 targets CD89 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FCAR Recombinant antibody, PBS Only (Capture) 84384-3-PBS (241426D1). 100 µg. Concentration 1 mg/mL.

Detection antibody: FCAR Recombinant antibody, PBS Only (Detector) 84384-2-PBS (241425F6). 100 µg. Concentration 1 mg/mL.

Alternative CD89 matched antibody pairs: MP01257-2, MP01257-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01257-1, CD89 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84384-3-PBS. Detection antibody: 84384-2-PBS. Standard: RP02338. Range: 0.156-10 ng/mL