

For Research Use Only

Mouse Ki67/MKI67 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01277-1

Capture Antibody Information

Catalog Number:
84432-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241768H7
Reactivity:
mouse
GenBank:
NM_001081117.2

Conjugate:
Unconjugated
Full name:
antigen identified by monoclonal antibody Ki 67
Gene ID:
17345

Detection Antibody Information

Catalog Number:
84432-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241768E2
Reactivity:
mouse
GenBank:
NM_001081117.2

Conjugate:
Unconjugated
Full name:
antigen identified by monoclonal antibody Ki 67
Gene ID:
17345

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01277-1 targets Ki67/MKI67 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Ki67/MKI67 Recombinant antibody, PBS Only (Capture) 84432-2-PBS (241768H7). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse Ki67/MKI67 Recombinant antibody, PBS Only (Detector) 84432-1-PBS (241768E2). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01277-1, MOUSE Ki67 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84432-2-PBS. Detection antibody: 84432-1-PBS. Standard: Eg1548. Range: 0.625-40 ng/mL.