

For Research Use Only

Mouse CEACAM-1/CD66a Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01297-1

Capture Antibody Information

Catalog Number:
84409-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241856E4
Reactivity:
mouse
GenBank:

Conjugate:
Unconjugated
Full name:
carcinoembryonic antigen-related
cell adhesion molecule 1
Gene ID:
26365

Detection Antibody Information

Catalog Number:
84409-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241856B1
Reactivity:
mouse
GenBank:

Conjugate:
Unconjugated
Full name:
carcinoembryonic antigen-related
cell adhesion molecule 1
Gene ID:
26365

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01297-1 targets CEACAM-1/CD66a in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CEACAM-1/CD66a Recombinant antibody, PBS Only (Capture/Detector) 84409-3-PBS (241856E4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Mouse CEACAM-1/CD66a Recombinant antibody, PBS Only (Detector) 84409-1-PBS (241856B1). 100 μ g. Concentration 1 mg/ml.

Alternative CEACAM-1/CD66a matched antibody pairs: MP00520-1, MP01297-2, MP01297-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01297-1, MOUSE CEACAM1/CD66a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84409-3-PBS. Detection antibody: 84409-1-PBS. Standard: Eg1423. Range: 0.313-40 ng/mL