

For Research Use Only

MLL4 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01298-1

Capture Antibody Information

Catalog Number:
84410-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241819A4
Reactivity:
human, rat
GenBank:
BC009337
Immunogen Catalog Number:
Ag35385

Conjugate:
Unconjugated
Full name:
myeloid/lymphoid or mixed-lineage leukemia 4
Gene ID:
9757

Detection Antibody Information

Catalog Number:
84410-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241819F4
Reactivity:
human
GenBank:
BC009337
Immunogen Catalog Number:
Ag35385

Conjugate:
Unconjugated
Full name:
myeloid/lymphoid or mixed-lineage leukemia 4
Gene ID:
9757

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA

Range:
0.781-100 ng/mL (Cytometric Bead Array)
46.9-3000 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01298-1 targets MLL4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array, Sandwich ELISA.

Capture antibody: MLL4 Recombinant antibody, PBS Only (Capture) 84410-1-PBS (241819A4). 100 µg. Concentration 1 mg/mL.

Detection antibody: MLL4 Recombinant antibody, PBS Only (Detector) 84410-2-PBS (241819F4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP01298-1, Human MLL4 Recombinant Matched Antibody Pair - PBS only. 84410-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag35385. 84410-2-PBS was HRP conjugated as the detection antibody. Range: 46.9-3000 pg/mL



The mean MLL4 concentration was determined to be 1.5 ng/mL in HeLa cell extract based on a 1.8 mg/mL extract load.



Cytometric bead array standard curve of MP01298-1, MLL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84410-1-PBS. Detection antibody: 84410-2-PBS. Standard: Ag35385. Range: 0.781-100 ng/mL