

ZNHIT2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01304-1

Capture Antibody Information

| | | |
|---|---|--|
| Catalog Number: 84418-3-PBS | Clone ID: 241868B8 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: zinc finger, HIT type 2 |
| Isotype: IgG | GenBank: BC052240 | Gene ID: 741 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag10351 | |

Detection Antibody Information

| | | |
|---|---|--|
| Catalog Number: 84418-4-PBS | Clone ID: 241868C11 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: zinc finger, HIT type 2 |
| Isotype: IgG | GenBank: BC052240 | Gene ID: 741 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag10351 | |

Applications

| | | |
|--|--|---|
| Tested Applications: Cytometric bead array | Range: 0.781-100 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|--|--|---|

Product Information

MP01304-1 targets ZNHIT2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ZNHIT2 Recombinant antibody, PBS Only (Capture) 84418-3-PBS (241868B8). 100 µg. Concentration 1 mg/ml.

Detection antibody: ZNHIT2 Recombinant antibody, PBS Only (Detector) 84418-4-PBS (241868C11). 100 µg. Concentration 1 mg/ml.

Alternative ZNHIT2 matched antibody pairs: MP01304-2, MP01304-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01304-1, ZNHIT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84418-3-PBS. Detection antibody: 84418-4-PBS. Standard: Ag10351. Range: 0.781-100 ng/mL.