

ZNHIT2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01304-2

Capture Antibody Information

Catalog Number: 84418-5-PBS	Clone ID: 241868E11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: zinc finger, HIT type 2
Isotype: IgG	GenBank: BC052240	Gene ID: 741
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10351	

Detection Antibody Information

Catalog Number: 84418-1-PBS	Clone ID: 241868A4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: zinc finger, HIT type 2
Isotype: IgG	GenBank: BC052240	Gene ID: 741
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10351	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01304-2 targets ZNHIT2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ZNHIT2 Recombinant antibody, PBS Only (Capture) 84418-5-PBS (241868E11). 100 µg. Concentration 1 mg/ml.

Detection antibody: ZNHIT2 Recombinant antibody, PBS Only (Detector) 84418-1-PBS (241868A4). 100 µg. Concentration 1 mg/ml.

Alternative ZNHIT2 matched antibody pairs: MP01304-1, MP01304-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage: Store at -80°C. The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer: PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01304-2, ZNHIT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84418-5-PBS. Detection antibody: 84418-1-PBS. Standard: Ag10351. Range: 0.781-100 ng/mL.