

For Research Use Only

CXCL9 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01308-2

Capture Antibody Information

Catalog Number:
84427-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241437G4

Reactivity:
human

GenBank:
BC095396

Immunogen Catalog Number:
Ag17932

Conjugate:
Unconjugated

Full name:
chemokine (C-X-C motif) ligand 9

Gene ID:
4283

Detection Antibody Information

Catalog Number:
84427-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241437E5

Reactivity:
human

GenBank:
BC095396

Immunogen Catalog Number:
Ag17932

Conjugate:
Unconjugated

Full name:
chemokine (C-X-C motif) ligand 9

Gene ID:
4283

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01308-2 targets CXCL9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CXCL9 Recombinant antibody, PBS Only (Capture) 84427-3-PBS (241437G4). 100 µg. Concentration 1 mg/ml.

Detection antibody: CXCL9 Recombinant antibody, PBS Only (Detector) 84427-2-PBS (241437E5). 100 µg. Concentration 1 mg/ml.

Alternative CXCL9 matched antibody pairs: MP01308-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01308-2, CXCL9/MIG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84427-3-PBS. Detection antibody: 84427-2-PBS. Standard: Ag17932. Range: 0.781-100 ng/mL