

NRP1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01310-2

Capture Antibody Information

Catalog Number:
84429-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241812A5

Reactivity:
human

GenBank:
BC007533

Conjugate:
Unconjugated

Full name:
neuropilin 1

Gene ID:
8829

Detection Antibody Information

Catalog Number:
84429-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241812C4

Reactivity:
human

GenBank:
BC007533

Conjugate:
Unconjugated

Full name:
neuropilin 1

Gene ID:
8829

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01310-2 targets NRP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NRP1 Recombinant antibody, PBS Only (Capture/Detector) 84429-1-PBS (241812A5). 100 µg. Concentration 1 mg/ml.

Detection antibody: NRP1 Recombinant antibody, PBS Only (Detector) 84429-3-PBS (241812C4). 100 µg. Concentration 1 mg/ml.

Alternative NRP1 matched antibody pairs: MP01310-1, MP01310-3, MP50479-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01310-2, Neuropilin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84429-1-PBS. Detection antibody: 84429-3-PBS. Standard: Eg2331. Range: 0.781-100 ng/mL.