

BAZ2A Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01313-1

Capture Antibody Information

Catalog Number:
82910-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230205H2
Reactivity:
human
GenBank:
NM_013449.4
Immunogen Catalog Number:
Ag33963

Conjugate:
Unconjugated
Full name:
bromodomain adjacent to zinc finger domain, 2A
Gene ID:
11176

Detection Antibody Information

Catalog Number:
82910-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230205H8
Reactivity:
human
GenBank:
NM_013449.4
Immunogen Catalog Number:
Ag33963

Conjugate:
Unconjugated
Full name:
bromodomain adjacent to zinc finger domain, 2A
Gene ID:
11176

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01313-1 targets BAZ2A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BAZ2A Recombinant antibody, PBS Only (Capture) 82910-5-PBS (230205H2). 100 µg. Concentration 1 mg/ml.

Detection antibody: BAZ2A Recombinant antibody, PBS Only (Detector) 82910-1-PBS (230205H8). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01313-1, BAZ2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82910-5-PBS. Detection antibody: 82910-1-PBS. Standard: Ag33963. Range: 0.625-20 ng/mL