

BAZ2A Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01313-2

Capture Antibody Information

Catalog Number:
82910-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230205F12
Reactivity:
human
GenBank:
NM_013449.4
Immunogen Catalog Number:
Ag33963

Conjugate:
Unconjugated
Full name:
bromodomain adjacent to zinc finger domain, 2A
Gene ID:
11176

Detection Antibody Information

Catalog Number:
82910-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230205F10
Reactivity:
human
GenBank:
NM_013449.4
Immunogen Catalog Number:
Ag33963

Conjugate:
Unconjugated
Full name:
bromodomain adjacent to zinc finger domain, 2A
Gene ID:
11176

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01313-2 targets BAZ2A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BAZ2A Recombinant antibody, PBS Only (Capture) 82910-4-PBS (230205F12). 100 µg. Concentration 1 mg/ml.

Detection antibody: BAZ2A Recombinant antibody, PBS Only (Detector) 82910-3-PBS (230205F10). 100 µg. Concentration 1 mg/ml.

Alternative BAZ2A matched antibody pairs: MP01313-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01313-2, BAZ2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82910-4-PBS. Detection antibody: 82910-3-PBS. Standard: Ag33963. Range: 1.25-20 ng/mL