

GINS2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01325-1

Capture Antibody Information

Catalog Number:
84486-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241870E9
Reactivity:
human
GenBank:
BC010164
Immunogen Catalog Number:
Ag9180

Conjugate:
Unconjugated
Full name:
GINS2 complex subunit 2 (Psf2 homolog)
Gene ID:
51659

Detection Antibody Information

Catalog Number:
84486-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241870B11
Reactivity:
human
GenBank:
BC010164
Immunogen Catalog Number:
Ag9180

Conjugate:
Unconjugated
Full name:
GINS2 complex subunit 2 (Psf2 homolog)
Gene ID:
51659

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01325-1 targets GINS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GINS2 Recombinant antibody, PBS Only (Capture) 84486-2-PBS (241870E9). 100 µg. Concentration 1 mg/ml.

Detection antibody: GINS2 Recombinant antibody, PBS Only (Detector) 84486-1-PBS (241870B11). 100 µg. Concentration 1 mg/ml.

Alternative GINS2 matched antibody pairs: MP01325-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01325-1, PSF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84486-2-PBS. Detection antibody: 84486-1-PBS. Standard: Ag9180. Range: 0.313-20 ng/mL