

For Research Use Only

Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01329-1

Capture Antibody Information

Catalog Number:
84460-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241725C1
Reactivity:
human
GenBank:
NM_001042724.2

Conjugate:
Unconjugated
Full name:
poliovirus receptor-related 2 (herpesvirus entry mediator B)
Gene ID:
5819

Detection Antibody Information

Catalog Number:
84460-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241725A10
Reactivity:
human
GenBank:
NM_001042724.2

Conjugate:
Unconjugated
Full name:
poliovirus receptor-related 2 (herpesvirus entry mediator B)
Gene ID:
5819

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01329-1 targets Nectin-2/CD112 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Nectin-2/CD112 Recombinant antibody, PBS Only (Capture) 84460-2-PBS (241725C1). 100 µg. Concentration 1 mg/ml.

Detection antibody: Nectin-2/CD112 Recombinant antibody, PBS Only (Detector) 84460-1-PBS (241725A10). 100 µg. Concentration 1 mg/ml.

Alternative Nectin-2/CD112 matched antibody pairs: MP01292-1, MP01292-2, MP01292-3, MP01329-2, MP01329-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01329-1, Nectin-2/CD112 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84460-2-PBS. Detection antibody: 84460-1-PBS. Standard: Eg1670. Range: 0.0781-50 ng/mL.