

TLR2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01332-1

Capture Antibody Information

Catalog Number:
84462-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241680E6

Reactivity:
human

GenBank:
BC033756

Conjugate:
Unconjugated

Full name:
toll-like receptor 2

Gene ID:
7097

Detection Antibody Information

Catalog Number:
84462-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241681G8

Reactivity:
human

GenBank:
BC033756

Conjugate:
Unconjugated

Full name:
toll-like receptor 2

Gene ID:
7097

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01332-1 targets TLR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TLR2 Recombinant antibody, PBS Only (Capture) 84462-1-PBS (241680E6). 100 µg. Concentration 1 mg/ml.

Detection antibody: TLR2 Recombinant antibody, PBS Only (Detector) 84462-4-PBS (241681G8). 100 µg. Concentration 1 mg/ml.

Alternative TLR2 matched antibody pairs: MP01332-2, MP01332-3, MP50388-1, MP50388-2, MP50388-3, MP50388-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01332-1, TLR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84462-1-PBS. Detection antibody: 84462-4-PBS. Standard: RP00009. Range: 0.781-100 ng/mL.