

For Research Use Only

ADAMTS2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01333-2

Capture Antibody Information

Catalog Number:
84476-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241860G5
Reactivity:
human
GenBank:
NM_014244
Immunogen Catalog Number:
Ag34957

Conjugate:
Unconjugated
Full name:
ADAM metalloproteinase with thrombospondin type 1 motif, 2
Gene ID:
9509

Detection Antibody Information

Catalog Number:
84476-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241860F12
Reactivity:
human
GenBank:
NM_014244
Immunogen Catalog Number:
Ag34957

Conjugate:
Unconjugated
Full name:
ADAM metalloproteinase with thrombospondin type 1 motif, 2
Gene ID:
9509

Applications

Tested Applications:
Cytometric bead array

Range:
0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01333-2 targets ADAMTS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ADAMTS2 Recombinant antibody, PBS Only (Capture/Detector) 84476-4-PBS (241860G5). 100 µg. Concentration 1 mg/mL.

Detection antibody: ADAMTS2 Recombinant antibody, PBS Only (Detector) 84476-3-PBS (241860F12). 100 µg. Concentration 1 mg/mL.

Alternative ADAMTS2 matched antibody pairs: MP01333-1, MP01333-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01333-2, ADAMTS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84476-4-PBS. Detection antibody: 84476-3-PBS. Standard: Ag34957. Range: 0.469-60 ng/mL.