

For Research Use Only



Mouse CD107b / LAMP2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01356-3

Capture Antibody Information

Catalog Number:
84474-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241893A7
Reactivity:
mouse
GenBank:
NM_001017959.2

Conjugate:
Unconjugated
Full name:
lysosomal-associated membrane
protein 2
Gene ID:
16784

Detection Antibody Information

Catalog Number:
84474-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241893E12
Reactivity:
mouse
GenBank:
NM_001017959.2

Conjugate:
Unconjugated
Full name:
lysosomal-associated membrane
protein 2
Gene ID:
16784

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP01356-3 targets CD107b / LAMP2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD107b / LAMP2 Recombinant antibody, PBS Only (Capture) 84474-5-PBS (241893A7). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Mouse CD107b / LAMP2 Recombinant antibody, PBS Only (Detector) 84474-2-PBS (241893E12). 100 μ g. Concentration 1 mg/ml.

Alternative CD107b / LAMP2 matched antibody pairs: MP01356-1, MP01356-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01356-3, MOUSE LAMP-2/CD107b Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84474-5-PBS. Detection antibody: 84474-2-PBS. Standard: Eg2029. Range: 0.625-80 ng/mL